


---

# Training for the adjustment of parabolic mirror at ELI-NP

---

Yoshihide Nakamiya  
LGED, ELI-NP/IFIN-HH



# Goals of the training course 1) Way of thinking for experimental research

---

## Basic cycle for experimental research

- Define your way scientifically (i.e. make hypothesis & reference)
- Execute your way
- Evaluate your way, that is, consider why your way is correct or wrong
- Modify your way if it is imperfect
- Conclude your study
- Outlook your study for the future

## Remarks

- We do not necessarily have the solution a priori or even the existence of the solution is not trivial for experimental research. So the basic cycle is important to tackle experimental research .
- The training coarse aims to experience the basic cycle in simple optical system.
- Keep “reasoning” when you take any actions.

# Goals of the training course 2) Contents

---

- Understand the basis of the alignment of optics and handling of optomechanics.
- Make your own configuration to align the optics on optical table.
- Do alignment of the mirrors, the beam expander, and focusing by off-axis parabolic mirror.
- Measure focal spot image by CMOS Camera
- Analyze the focal spot image.
- Present your configuration, the way of methods, and your results

# Schedule of training course

---

- Kick-off meeting for introduction of training course **(Today)**
  - Introduction of training course (basis of alignment, goals and laboratory works)
  - Explanation of equipment at LP14
- 5 days for laboratory works for each group (in three groups) at lp14
  - Start-up discussion **at 10:00 on Monday**
  - Finish lab. works **by 15:00 on Friday**
- Presentation by each group in **30 July (Wed.)**
- Back-up day in **31 July (Thu.)**

Day	Time	Contents	Remarks
4 July (Fri)	09:30 - 12:00	Introduction of training course	
	13:00 - 15:00	Introduction of laboratory	Explanation of equipment/optomechanics/tools at LP14
07 July - 11 July	-	Training at LP14 for Group A	Finish by 15:00 on the final day (Friday) for next group
14 July - 18 July	-	Training at LP14 for Group B	
21 July - 25 July	-	Training at LP14 for Group C	
30 July (Wed)	09:30 - 11:30	Presentation by Group A and B	30 minutes presentation + 20 minutes discussion
	13:00 - 14:30	Presentation by Group C	
	14:30 - 17:00	Summary and dicussion	
31 July (Thu)	09:30 - 12:00	Back-up day for further discussion	
4 Aug - 8 Aug	-	Back-up week	

# Schedule of group activity in laboratory

Group	Period of Lab activities	Name	Department
A	07 July (Mon) - 11 July (Fri)	Annemarie Cristmann	LDED
		Fattima	LSD
		Popescu Vlad	LSD
		Saidbek	LSD
B	14 July (Mon) - 18 July (Fri)	Cana Alexandra Cristina	GSD
		Claudia Jerca	LDED
		Diana Naum	LDED
		Iancu Andrei Daniel	GSD
		Ioan Paul Pârlea	GDED
		Ionut Slabu	LGED
C	21 July (Mon) - 25 July (Fri)	Cristi Alexe	LSD
		Diana Catana	LDED
		Emilia Matache	LSD
		Grigore Alexandra-Elena	LSD
		Octavian Cretu	LDED
		Patrascoiu Robert	GSD
		Sandu Denis	LSD


# Remarks for 1 week laboratory works by each group

---

- Start-up discussion **at 10:00 on Monday**
  - Brief remarks of equipment/tools/optomechanics
  - Confirmation of safety issues
- Completion of laboratory work **at 15:00 on Friday**
  - Deconstruct the system for initialization for next group after 15:00
- You should concentrate on taking the data and collecting the information of your system on Friday, regardless of your progress.

# Block diagram of the optical system in the lab.

---



- All the optical components are installed on the mounts and mechanics.
- The training coarse starts from to consider your optical setup and align them almost from the scratch.
  - We skipped the process of considering/designing/finding each components due to the limitation of sources. Just keep in your mind when you execute your experimental researches

# Definition of terminology to adjust mirror mount





illustration is cited from Wikipedia  
[https://en.wikipedia.org/wiki/Aircraft\\_principal\\_axes#/media/File:Yaw\\_Axis\\_Corrected.svg](https://en.wikipedia.org/wiki/Aircraft_principal_axes#/media/File:Yaw_Axis_Corrected.svg)

“Yaw” is rotation horizontally  
“Pitch” is rotation vertically



Typical kinematic mirror mount


Knob for Pitch

Knob for Yaw




# Bais of alignment

## Principle



- (0) Initial ray vector
- (1) “Offset” alignment  
(equiv. to parallel shift)
- (2) “Pointing” alignment  
(equiv. to angle correction)


## Alignment



- If M2 has adjustability of
- horizontal & vertical shift (by translation stage)
  - yaw & pitch (by mirror mount)
- + if  $L1=0$ ,  $L2=\infty$ ,  
M2 adjustment works exactly as principle


# Alignment of two mirrors

---



- Yaw and Pitch for M1 is used for “offset” alignment with I1.
- Yaw and Pitch for M2 is used for “pointing (angle)” alignment with I2
- Since  $L1 \neq 0$  and  $L2 \neq \infty$ , we need iteration of alignment with finite precision of angle
- Since Adjustment of M1 and M2 are not independent operation, we need iteration of alignment

# Alignment of beam expander



- Beam expander is made of a Keplerian telescope configuration by two convex lens.
- Install CV1 and CV2 and optimize the lens position to keep reference laser alignment
- Optimize laser divergence by adjusting the distance between CV1 and CV2
- Collimated beam is necessary as output

# Property of paraboloid

Parabolic Mirror





Illustration is cited from Thorlabs page:  
[https://www.thorlabs.co.jp/newgrouppage9.cfm?objectgroup\\_id=14211](https://www.thorlabs.co.jp/newgrouppage9.cfm?objectgroup_id=14211)

- Key property for focusing by parabolic mirror
- Collimated incoming laser (well control of divergence)
  - Sensitive to angle of incident of incoming laser (good precision of angle of incidence)
  - Insensitive of offset shift of incoming laser


# Off-axis parabolic mirror

Parabolic Mirror




- If one cut paraboloid according to light blue line  
-> On-axis parabolic mirror
- If one cut paraboloid according to light green line  
-> Off-axis parabolic mirror
- Angle of Off-axis parabolic mirror is defined as  $\theta$  (Full-angle) or  $\theta/2$  (AOI, angle of incidence) in specification

Parabolic Mirror




# Coarse alignment of off-axis parabolic mirror (OAP)



- Alignment of incident angle to OAP
- Alignment of outgoing angle from OAP to set initial position to bring the laser to expected focal-spot point
- Install CMOS camera + ND5 (See photograph)
- Check the astigmatism on image

# Tips: Adjustability of off-axis parabolic mirror and CMOS Camera

---



## 5 axes for OAP

- V (Vertical translation)
- H (Horizontal translation)
- L (Longitudinal translation)
- Pitch
- Yaw

## 3 axes for CMOS camera

- V (Vertical translation)
- H (Horizontal translation)
- L (Longitudinal translation)





# Astigmatism by off-axis parabolic mirror





Image at A



Image at O






- Not focused well
- Non-trivial shape  
(white circle show expected optimal focusing image)

Image at B






# Astigmatism by off-axis parabolic mirror

Astigmatic case (misaligned at a few degrees)




As same as (e),  
But gain is optimal  
(no saturation)

Little astigmatic case




# Confirmation of “input” laser quality if your beam is astigmatic



Convex Lens

Aspheric objective lens



- You can check input laser beam quality by focusing lens
- You may need aspheric lens to see detailed shape for focusing due to necessity of magnification




# Fine adjustment of OAP for shaping focal-spot profile

---




- Install aspheric lens on CMOS camera to find focal point
- Adjust Yaw & H to keep the image on the camera center
  - To get the circular shape
  - To get the smallest spot size
  - To get the highest intensity spot
- (If necessarily) Adjust Pitch & V to keep the image on the camera center to get the smallest spot size as well as the highest intensity spot.

# Calibration of magnification



Move the camera by knob and see how much intensity centroid (pixels) moves in on sensor plane

# Focal-spot measurement and evaluation



## Spot radius $w_0$ of Gaussian laser beam

$$w_0 \sim \frac{2}{\pi} \lambda f_{num},$$

$$f_{num} = \frac{f}{\Phi},$$

$\lambda$ : wavelength

$f_{num}$  : f number

$f$ : focal distance

$\Phi$ : laser-beam diameter before focusing

# Missions of your laboratory activities

---

- Alignment of mirrors and beam expander
- Control of divergence of beam expander (Collimated laser beam)
- Alignment towards OAP (angle of incidence to OAP)
- Coarse alignment of OAP (angle of outgoing from OAP)
- Measure focal image without aspheric objective lens (astigmatism check)
- Alignment of CMOS camera with aspheric objective lens
- Fine alignment of OAP (optimal shaping of focal-spot profile)
- Data acquisition 1) Focal-image by focusing lens before OAP (with calibration)
- Data acquisition 2) Focal-image by OAP after OAP (with calibration)
- Offline data analysis
  - Focal-spot size with physical scale ( $\mu\text{m}$ )
  - Evaluate how good or bad your focusing image is
- Make your presentation
- (Advanced) play further to see tendency (e.g. detune worse direction on purpose to see how astigmatism grows etc.)



# Presentation by each group in 30 July (Wed.)

---

The meeting for presentation of your laboratory activities.

- Presentation in **30 July (Wed.)**
- Each group should **present for 30 minutes + discussion for 20 minutes**,

Your presentation should include the following points

- Show your optical setup
- Show your way of the alignment. In particular,
  - Optimization of beam expander
  - Optimization of angle of incident to OAP
  - Optimization of angle of outgoing from OAP
- Show your measurement, analysis and result
  - Data sets
  - The way of analysis
  - Analysis results (Calibration and spot size)
- Show your opinions/insights about, for example,
  - What is the most difficult/critical part for your activities ?
  - How do you tackle and take care of the difficulties ?
  - Any strange/interesting phenomenon you noticed
- Show your Summary
  - What is your achievement ?
  - What is remaining issues and improvement ?
  - Anything else ?



# Laser safety

---

- Never do your eye level to be the same on the laser beam
- Remove the reflective accessory (for instance, bracelet, watch etc.)
- Build the beam blocker at the termination in your alignment process, in particular, for the direction to the area where the other people access.
- Wear the safety goggles

# Short discussions after the meeting

---

Confirm grouping and group members

Confirm equipment/optomechanics/tools

# Back up

---

# Component List

Name	Qty	Model	Manufacturer	Remarks
HeNe Laser	1	HNL050R	Thorlabs	
Adjustable ND filter	1	NDC-50C-4M-B	Thorlabs	
1 inch Ag mirror	3	PF10-03-P01	Thorlabs	
1 inch Al mirror	1	PF10-03-G01	Thorlabs	
1 inch Al mirror	1	092-3015	Eksma	
2 inch Ag mirror	1	PF20-03-P01	Thorlabs	
2 inch Au mirror	1	20Z40ER.4	Newport	
20 mm Iris	1	IDA20/M	Thorlabs	Minimal aperture size is 1.2 mm
25 mm Iris	1	ID25SS/M	Thorlabs	Minimal aperture size is 1.4 mm
37 mm Iris	1	ID37Z/M	Thorlabs	Minimal aperture size is 0 mm
50 mm Iris	1	ID50Z/M	Thorlabs	Minimal aperture size is 0 mm
Convex f10	1	AC080-010-A-ML	Thorlabs	Achromatic
Convex f100	1	SLB-50.8-100PM	Thorlabs	Plano-convex
Convex f150	1	SLB-50.8-150PM	Thorlabs	Plano-convex
Off-axis parabolic	1	35516	Edmund	Off-set full angle 60 deg., EFL=101.6mm
Aspheric lens f4.5 mm	1	C230TMD-B	Thorlabs	Mount on CMOC camera
1 inch ND	1	NE50B	Thorlabs	Mount on CMOC camera
CMOS camera	1	acA2440-20gm	Basler	